

## Annex 1

# Instruction for Use

## THYROXINE (T4) ELISA KIT

Catalog No. E1002

### **INTENDED USE**

The Autobio T4 assay is designed for the quantitative determination of Thyroxine (T4) concentration in human serum.

### **INTRODUCTION**

Thyroxine or 3, 5, 3', 5'-tetraiodo-L-thyronine (T4) is the major hormone produced by the thyroid gland. It has a molecular weight of 777 daltons and is synthesized by iodination of tyrosine residues on thyroglobulin. Proteolytic cleavage of follicular thyroglobulin releases T4 into the bloodstream. Greater than 99% of T4 is reversibly bound to three plasma proteins in blood - thyroxine binding globulin (TBG) binds 70%, thyroxine binding pre-albumin (TBPA) binds 20%, and albumin binds 10%. Approximately 0.03% of T4 is in the free, unbound state in blood at any time.<sup>1,2</sup>

Diseases effecting thyroid function may present a wide array of confusing symptoms. Measurement of total T4 by immunoassay is the most reliable and convenient screening test available to determine the presence of thyroid disorders in patients. Increased levels of T4 have been found in hyper-thyroidism due to Grave's disease and Plummer's disease and in acute and subacute thyroiditis. Low levels of T4 have been associated with congenital hypothyroidism, myxedema, chronic thyroiditis (Hashimoto's disease), and with some genetic abnormalities.<sup>3,4,5</sup>

### **PRINCIPLE OF THE TEST**

In the T4 EIA, a certain amount of anti-T4 antibody is coated on microtiter wells. A measured amount of patient serum, and a constant amount of T4 conjugated with horseradish peroxidase are added to the microtiter wells. 8-anilino-1-naphthalene sulfonate (ANS) is used to displace T4 from proteins to enable the measurement of total circulating T4. During incubation, T4 and conjugated T4 compete for the limited binding sites on the anti-T4 antibody. After 60 minutes incubation at 37°C, the wells are washed by Wash Solution. Then Substrate Solution and Chromogen Solution are added and incubated for 20 minutes, resulting in the development of blue color. The color development is stopped with the addition of Stop Solution, and the color is changed to yellow and absorbance is measured spectrophotometrically at 450 nm. The color intensity is inversely related to the concentration of T4 in the test sample.

### **MATERIALS PROVIDED**

1. Antibody Coated Microtiter Plate: Microplate with anti-T4 Antibody coated wells (1 plate, 96 wells)
2. Enzyme Conjugate Concentrate: Horseradish Peroxidase (HRP) labeled T4 in Stabilizing Buffer (1 vial, 1.2 ml, 20X)
3. Enzyme Conjugate Diluent ( 2 vials, 7.5 ml/ea)
4. Reference Standards: 0, 1.0, 2.5, 5.0, 15, 30µg/dl T4 in HEPES solution with preservatives. 1µg/dl of the reference standard is equivalent to 1µg/dl standard of the National Institute for the Control of Pharmaceutical and Biological Products (NICPBP), P. R.C. (6 vials, 1ml/ea)
5. Wash Solution Concentrate: PBS-Tween (1 bottle, 25 ml, 40X)
6. Substrate Solution: Hydrogen Peroxide (1 vial, 7.5 ml)
7. Chromogen Solution: Tetramethylbenzidine (1 vial, 7.5 ml)
8. Stop Solution: 1.0M H<sub>2</sub>SO<sub>4</sub> (1 vial, 7.5 ml)

### **MATERIALS NOT PROVIDED**

The following materials are required but not provided in the kit:

1. Distilled water
2. Precision pipettes for delivery of 20-200µl, 100-1000µl (the use of accurate pipettes with disposable plastic tips is recommended)
3. Microplate ELISA reader with a bandwidth of 10nm or less and an optical density range of 0-3.5 OD or greater at 450nm wavelength

4. Magnetic stirrer
5. Vortex Mixer or equivalent
6. Washer for microplates
7. Quality control specimens
8. Incubator
9. Absorbent paper

#### **STORAGE OF TEST KIT AND INSTRUMENTATION**

1. Unopened test kits should be stored at 2-8°C upon receipt and the microplate should be kept in a sealed bag with desiccants to minimize exposure to damp air. The test kit may be used throughout the expiration date of the kit (One year from the date of manufacture). Refer to the package label for the expiration date.
2. Opened test kits will remain stable for at least two months, provided it is stored as prescribed above.

#### **SPECIMEN COLLECTION AND PREPARATION**

1. Serum is the recommended sample type for this assay. Plasma samples collected in tubes containing EDTA, heparin, or oxalate may interfere with test procedures and should be avoided.
2. Collect all blood samples observing universal precautions for venipuncture.
3. Allow samples to clot for 1 hour before centrifugation.
4. Avoid grossly hemolytic, lipemic or turbid samples.
5. Prior to use, specimens should be capped and stored up to 48 hours at 2-8°C. For longer storage, freeze the specimens at -20°C. Thawed samples must be mixed prior to testing.

#### **PRECAUTIONS AND WARNINGS**

1. For *in vitro* diagnostic use only.
2. Handling of reagents, serum specimens should be in accordance with local safety procedures.
3. All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.
4. Avoid any skin contact with all reagents, Stop Solution contains H<sub>2</sub>SO<sub>4</sub>, in case of contact, wash thoroughly with water.
5. Sodium azide in reference standards can react with copper and lead plumbing to form explosive metal azides. On disposal, flush reagents with a large volume of water to prevent the buildup of azides, if disposal into a drain is in compliance with federal, state, and local requirements.
6. Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

#### **REAGENT PREPARATION**

1. All reagents should be allowed to reach room temperature (18-25 °C) before use.
2. Adjust the incubator to 37°C.
3. Prepare Enzyme Conjugate Reagent: add 0.1 ml of Enzyme Conjugate Concentrate to 2.0 ml of Enzyme Conjugate Diluents (1:20 dilution), and mix well. The amount of conjugate diluted is depending on your assay size. The Conjugate Reagent is stable at 4°C for 7 days.
4. Prepare Wash Solution: add 25ml of Wash Solution Concentrate to 1000ml of distilled water, and mix well with magnetic stirrer. The Wash Solution is stable at R.T. for two months.

#### **IMPORTANT NOTES**

1. Do not use reagents after expiration date.
2. Do not mix or use components from kits with different lot numbers.
3. It is recommended that no more than 32 wells be used for each assay run, if manual pipette is used, since pipetting of all standards, specimens and controls should be completed within 5 minutes. A full plate of 96 wells may be used if automated pipette is available.
4. Replace caps on reagents immediately. Do not switch caps.
5. The wash procedure is critical. Insufficient washing will result in poor precision and falsely elevated absorbance readings.
6. When the Enzyme Conjugate Concentrate is diluted in the Enzyme Conjugate Diluent, the obtained Conjugate Reagent is stable at 4°C for 7 days. We recommend the amount of Enzyme Conjugate Concentrate diluted is depending on your assay size.

**ASSAY PROCEDURE**

1. Secure the desired number of coated well in the holder. Make data sheet with sample identification.
2. Dispense 50 µl of standards, samples, and controls into appropriate wells.
3. Dispense 100 µl of Enzyme Conjugate Reagent into each well.
4. Thoroughly mix for 30 seconds. It is important to have complete mixing in this step.
5. Incubate at 37°C for 60 minutes.
6. Remove the incubation mixture by flicking plate contents into a waste container.
7. Rinse and flick the micro plate 5 times with wash solution.
8. Strike the wells sharply onto absorbent paper to remove residual water droplets.
9. Dispense 50 µl Chromogen Solution into each well.
10. Dispense 50 µl Substrate Solution into each well. Gently mix for 15 seconds.
11. Incubate at room temperature in the dark for 20 minutes without shaking.
12. Stop the reaction by adding 50 µl of Stop Solution to each well.
13. Gently mix for 15 seconds. It is very important to make sure that the blue color changes to yellow color completely.
14. Read absorbance at 450nm with a microplate reader within 15 minutes.


**CALCULATION OF RESULTS**

1. Calculate the mean value from any duplicate reagents. Where appropriate, the mean values should be used for plotting.
2. On linear graph paper plot the absorbance values (ordinate) for each calibrator against the corresponding concentration of T4 in µg/dl (abscissa) and draw a calibration curve through the calibrator points by connecting the plotted points with straight lines.
3. Read the concentration for each control and sample by interpolation on the calibration curve.
4. Computer assisted data reduction will simplify these calculations. If automatic result processing is used, a 4-parameter logistic function curve fitting is recommended.

**EXAMPLE OF STANDARD CURVE**

A typical standard curve shown below is for the purpose of illustration only, and should never be used instead of the real time calibration curve.

| T4 (µg/dl) | Absorbance (450nm) |
|------------|--------------------|
| 0          | 3.015              |
| 1.0        | 1.852              |
| 2.5        | 1.115              |
| 5.0        | 0.652              |
| 15         | 0.300              |
| 30         | 0.176              |



**EXPECTED VALUES**

Each laboratory should establish its own normal range. These values are given only for guidance.

|                       |      |
|-----------------------|------|
| Sample Number         | 200  |
| Average Value (µg/dl) | 8.95 |



|   |          |
|---|----------|
| Standard Deviation ( $\sigma$ )                   | 2.01     |
| Normal Range ( $\pm 2\sigma$ , $\mu\text{g/dl}$ ) | 5.0-13.0 |

**PERFORMANCE**

**A. Sensitivity**

Twenty zero standards were assayed along with a set to other standards. The detection limit, defined as the apparent concentration corresponding to two standard deviations below the average absorbance at zero binding, was not higher than 0.4 $\mu\text{g/dl}$ .

**B. Specificity**

The cross-reactivity of the total T4 assay kit with T3 and rT3 was determined by adding these hormones to zero standards. The color intensity produced was then determined.

| Interferent | Concentration | Measured Value ( $\mu\text{g/dl}$ ) | Crosstalk Rate (%) |
|-------------|---------------|-------------------------------------|--------------------|
| T3          | 500ng/ml      | 1.06                                | 2.12               |
| rT3         | 500ng/ml      | 0.84                                | 1.68               |

**C. Precision**

**a. Intra-Assay Precision**

Intra-Assay Precision was determined by assaying 20 replicates of each of 2 sera; low and high.

| Serum | Number | Mean  | SD   | RSD (%) |
|-------|--------|-------|------|---------|
| Low   | 20     | 3.92  | 0.22 | 5.65    |
| High  | 20     | 11.12 | 0.40 | 3.58    |

**b. Inter-Assay Precision**

Inter assay Precision was determined by assaying duplicates of 2 serum pools in 20 separate runs, using a standard curve constructed for each run.

| Serum | Number | Mean  | SD   | RSD (%) |
|-------|--------|-------|------|---------|
| Low   | 20     | 3.84  | 0.37 | 9.64    |
| High  | 20     | 10.74 | 1.08 | 10.06   |

**D. Accuracy**

For 90 samples in the range of 1.5 $\mu\text{g/dl}$  to 25 $\mu\text{g/dl}$ , the relationship between the Autobio T4 ELISA Test and the Biocheck T4 ELISA Test is described by the equation:

| Reference        | No. of Specimens | Least Square Regression Analysis | Correlation Coefficient |
|------------------|------------------|----------------------------------|-------------------------|
| Biocheck (ELISA) | 90               | $Y=1.0776x + 1.2921$             | 0.954                   |

**LIMITATIONS**

Reliable and reproducible results will be obtained when the assay procedure is carried out with a complete understanding of the package insert instructions and with adherence to good laboratory practice.

Serum T4 concentration is dependent upon a multiplicity of factors: hypothalamus gland function and its regulation, thyroxine binding globulin (TBG) concentration, and the binding of T4 to TBG. Thus, total T4 concentration alone is not sufficient to assess clinical status.









Heterophilic antibodies in human serum can react with reagent immunoglobulin, interfering with *in vitro* immunoassays. Patients routinely exposed to animals or to animal serum products can be prone to this interference thus anomalous values may be observed. Additional information may be required for diagnosis.

For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

**QUALITY CONTROL**

Good laboratory practice requires that quality control specimens be run with each calibration curve to verify assay performance. To assure proper performance, a statistically significant number of controls should be assayed to establish mean values and acceptable ranges. Controls containing sodium azide should not be used.

**SYMBOLS**

|   |  |
|---|--|
|    | <b>BATCH CODE</b>                              |
|    | <b>USE BY</b>                                  |
|    | <b>MANUFACTURER</b>                            |
|    | <b>CONTAINS SUFFICIENT FOR &lt;n&gt; TESTS</b> |
|    | <b>IN VITRO DIAGNOSTIC MEDICAL DEVICE</b>      |
|    | <b>TEMPERATURE LIMITATION</b>                  |
|   | <b>CATALOGUE NUMBER</b>                        |
|  | <b>CONSULT INSTRUCTIONS FOR USE</b>            |

**REFERENCES:**

- [1]Young DS, Pestaner LC, and Gilberman. Effects of Drugs on Clinical Laboratory Test. Clinical Chemistry. 21, 3660 (1975)
- [2]Sterling L, Diagnosis and treatment of Thyroid Disease, Cleveland, CRC Press, 9, 51 (1975)
- [3]Surk M L, Chopra I J, Mariash C N. et al. American thyroid association guidelines for use of laboratory tests in thyroid disorders [J]. JAMA, 1990, 263: 1529—1532.
- [4]Glinoe D. The regulation of thyroid function in pregnancy: path ways of endocrine adaptation from physiology to pathology [J]. Endocrinol Rev. 1997. 18: 404—33.
- [5]Demem L M, Spencer C A. Laboratory medicine practice guidelines for laboratory support for the diagnosis and monitoring of thyroid diseases [J]. National Academy of Clinical Biochemistry, 2002.

For order and inquires, please contact



AUTOBIO DIAGNOSTICS CO., LTD.  
 ADD: No.87 Jingbei Yi Road, National Eco & Tech Development Area,  
 Zhengzhou , China 450016  
 Tel: +86-371-67985313 Fax: +86-371-67985804  
 Web: www.autobio.com.cn