

---

## INSULIN (INS) CHEMILUMINESCENCE IMMUNOASSAY KIT

Catalog No. CL1201-2

### **INTENDED USE**

The Autobio insulin (INS) chemiluminescence immunoassay (CLIA) kit is intended for the quantitative determination of INS concentration in human serum.

### **INTRODUCTION**

Insulin is a member of structurally related regulatory proteins; other proteins in this group include the insulin-like growth factors and relaxin. Insulin is produced by the  $\beta$ -cells of the pancreatic islets and is initially synthesized as a 12kDa preprohormone, which undergoes intracellular processing to a 9kDa, 86-amino acid prohormone and subsequent packaging in storage granules. Within these granules, disulfide bonds are formed between the A and B chains of the insulin molecule and the C-peptide region is cleaved, resulting in the 51-amino acid, 6kDa mature insulin molecule. Upon stimulation, the islet cells release equimolar amounts of insulin and C-peptide, and small amounts of proinsulin and other intermediates (<5% of normal total insulin secretion).

Basal and glucose-stimulated circulating insulin concentrations are relatively stable during infancy and childhood, and increase during puberty due to decreased insulin sensitivity. Insulin concentrations tend to be higher in obese individuals, particularly those with an increased proportion of visceral (abdominal) fat. Glucose counter-regulatory hormones, such as glucagons, gluconcorticoids, growth hormone and epinephrine, decrease insulin sensitivity and action; insulin levels may increase during exogenous administration of these substances.

Measurement of circulating insulin concentrations may be useful in the diagnostic evaluation of several conditions. Elevated serum insulin levels in the presence of low glucose concentrations may be indicative of pathologic hyperinsulinism, e.g. nesidioblastosis and islet-cell tumor. Elevated serum fasting insulin levels with normal or elevated glucose concentrations, and exaggerated insulin and glucose response to exogenous glucose administration are characteristic of the insulin-resistant forms of glucose intolerance and diabetes mellitus and other insulin resistant conditions. High circulating insulin concentrations may be involved in the pathogenesis of hypertension and cardiovascular disease. Conversely, low insulin concentrations in the presence of hyperglycemia suggest insulin-deficiency, e.g. insulin-dependent or Type I diabetes mellitus. Measurement of immediate or first-phase insulin secretion after an acute glucose load may be predictive of Type I diabetes mellitus.

Endogenous anti-insulin antibodies may be observed in the prediabetic and symptomatic phases of insulin-dependent diabetes mellitus, presumably due to an autoimmune disorder. Anti-insulin antibodies are also commonly seen in insulin-treated individuals such antibodies can interfere with immunoassays for insulin, and techniques have been described to remove the antibodies prior to assay for "free" insulin. The presence of excess proinsulin and genetic insulin variants can also interfere with insulin estimations.

### **PRINCIPLE OF THE TEST**

The INS CLIA test is a solid phase two-site immunoassay. One monoclonal antibody is coated on the surface of the microtiter wells and another monoclonal antibody labeled with horseradish peroxidase is used as the tracer. The INS molecules present in the standard solution or serum are "sandwiched" between the two antibodies. Following the formation of the coated antibody-antigen-antibody-enzyme complex, the unbound antibody enzyme labels are removed by washing. The horseradish peroxidase activity bound in the wells is then assayed by chemiluminescence reactions. The Related Light Unit (RLU) of the reaction is proportional to the concentration of INS present in the sample.

### **MATERIALS PROVIDED**

1. Antibody Coated Microtiter Plate: Microplate with monoclonal antibodies to insulin (anti-INS MAb) coated wells (1 plate, 48wells/96 wells)
2. Enzyme Conjugate Reagent: Horseradish peroxidase (HRP) labeled anti-INS (MAb) in Stabilizing Buffer (1 vial, 3ml/6.0ml)
3. Reference Standards: 5, 20, 50, 100, 200 IU/ml (5 vials, Lyophilized)
4. Substrate A: (1 vial, 3.5ml/6.0ml)

5. Substrate B: (1 vial, 3.5ml/6.0ml)
6. PBS-T Powder: PBS-Tween (1 bag, 5g)

**MATERIALS NOT PROVIDED**

The following materials are required but not provided in the kit

1. Precision pipettes and tips, 0.05ml, 0.1ml, 1.0ml
2. Distilled water
3. Vortex mixer
4. Absorbent paper or paper towel
5. Graph paper
6. Luminometer

**STORAGE OF TEST KIT AND INSTRUMENTATION**

1. Unopened test kits should be stored at 2~ 8°C upon receipt and the microtiter plate should be kept in a sealed bag with desiccants to minimize exposure to damp air. The test kit may be used throughout the expiration date of the kit. Refer to the package label for the expiration date.
2. Opened test kits will remain stable until the expiring date shown, provided it is stored as prescribed above.

**SPECIMEN COLLECTION AND PREPARATION**

1. Blood should be drawn using standard venipuncture techniques and the serum should be separated from the red blood cells as soon as practical. Avoid grossly hemolytic, lipemic or turbid samples.
2. Plasma samples collected in tubes containing EDTA, heparin, or oxalate may interfere with test procedures and should be avoided.
3. Specimens should be capped and may be stored up to 48 hours at 2 ~ 8°C, prior to assaying. Specimens held for a longer time can be frozen at -20°C. Thawed samples must be mixed prior to testing.

**PRECAUTIONS AND WARNINGS**

1. For *in vitro* diagnostic use only.
2. Handling of reagents, serum specimens should be in accordance with local safety procedures.
3. The reference standards contain human source components, which have been tested and found non-reactive for hepatitis B surface antigen as well as HIV antibody. All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, the standards and components containing animal substances should be treated as potentially infectious.
4. Avoid any skin contact with all reagents.
5. Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

**REAGENT PREPARATION**

1. All reagents should be brought to room temperature (18 ~ 25°C) before use. All reagents should be mixed through gently inverting or swirling prior to use. Do not induce foaming.
2. To prepare Washing Buffer: add 1 bag of washing buffer concentration to 500ml of distilled water, and mix well. The washing buffer is stable at room temperature at least for two weeks.

**IMPORTANT NOTES**

1. The wash procedure is critical. Insufficient washing will result in poor precision and falsely elevated RLU values.
2. It is recommended that no more than 32 wells be used for each assay run, if manual pipetting is used, since pipetting of all standards, specimens and controls should be completed within 5 minutes. A full plate of 96 wells may be used if automated pipetting is available.
3. Duplication of all standards and specimens, although not required, is recommended.

**ASSAY PROCEDURE**

1. Secure the desired number of coated wells in the holder. Dispense 50µl of reference standards, specimens, and controls into appropriate wells.
2. Dispense 50µl of Enzyme Conjugate Reagent to each well. Mix gently for 30 seconds.
3. Incubate at 37°C for 60 minutes.

4. At the end of the incubation, remove the incubation mixture by emptying the plate content into a waste container. Rinse and empty the microtiter plate 5 times with washing buffer. Strike the microtiter plate sharply onto absorbent paper or paper towels to remove all residual water droplets. The volume of each well is about 350µl.
5. Dispense 50µl of substrate A, then 50µl of substrate B into each well. Gently mix for 10 seconds.
6. Incubate at room temperature in the dark for 5 minutes without shaking and read the RLU values with a Luminometer.


**CALCULATION OF RESULTS**

1. Calculate the mean value from any duplicate reagents. Where appropriate, the mean values should be used for plotting.
2. Plot the log<sub>10</sub>RLU for each reference standard against the common logarithm of corresponding concentration of INS in mIU/ml on logarithmic graph paper, with RLU values on the Y-axis and concentration on the X-axis.
3. Read the concentration for each control and sample by interpolation on the calibration curve.
4. Computer assisted data reduction will simplify these calculations. If automatic result processing is used, a linear regression logistic function curve fitting is recommended.

**EXAMPLE OF STANDARD CURVE**

A typical standard curve shown below is for the purpose of illustration only, and should never be used instead of the real time calibration curve.

| INS ( IU/ml) | RLU     |
|--------------|---------|
| 5            | 561.2   |
| 20           | 2885.2  |
| 50           | 8186.5  |
| 100          | 18935.5 |
| 200          | 41102.3 |



**EXPECTED VALUES**

Each laboratory must establish its own normal ranges based on patient population. The normal range is between 2.0 IU/ml and 20 IU/ml, which were determined by testing 200 fasting sera samples.

**PERFORMANCE**

**A. Sensitivity**

The lower detection limit is calculated from the standard curve by identifying the concentration corresponding to the mean RLU of standard diluent (based on 10 replicate analyses) plus 2 SD. Therefore, the sensitivity of the Autobio INS CLIA kit is estimated to be not higher than 2.5 IU/ml.

**B. Specificity**

No interference was detected with the performance of Autobio INS CLIA upon addition of massive amounts of the following substances to a human serum pool.

| Interferen\$ | Concentration |
|--------------|---------------|
| Proinsulin   | 1000pmol/L    |

**C. Precision**

**a. Intra-assay Precision**

Intra-assay Precision was determined by assaying 20 replicates of each control sera.

| Serum      | Number | Mean  | SD   | CV (%) |
|------------|--------|-------|------|--------|
| Low titer  | 20     | 16.26 | 0.83 | 5.10   |
| High titer | 20     | 86.58 | 2.03 | 2.34   |

**b. Inter-assay Precision**

Inter-assay precision was determined by assaying duplicates of each control sera in 10 separate runs.

| Serum      | Number | Mean  | SD   | CV (%) |
|------------|--------|-------|------|--------|
| Low titer  | 10     | 15.38 | 1.10 | 7.12   |
| High titer | 10     | 81.50 | 2.73 | 3.36   |

**D. High Dose Hook Effect**

In this INS assay, patient samples spiked to INS levels up to 1000 IU/ml do not demonstrate a paradoxical decrease in the RLU's (high dose hook effect).

**E. Accuracy**

For samples in the range of 2.5 IU/ml to 200 IU/ml, the correlation between the Autobio INS CLIA kit and Beckman access 2<sup>®</sup> INS assay, are described by the equation:

| Reference                     | Number of Specimens | Least Square Regression Analysis | Correlation Coefficient |
|-------------------------------|---------------------|----------------------------------|-------------------------|
| Beckman access 2 <sup>®</sup> | 180                 | $y = 0.9438x + 2.0451$           | 0.9715                  |









**LIMITATIONS**

1. Reliable and reproducible results will be obtained when the assay procedure is carried out with a complete understanding of the package insert and with adherence to good laboratory practice.
2. Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays. Patients routinely exposed to animals or to animal serum products can be prone to this interference thus anomalous values may be observed. Additional information may be required for diagnosis.
3. Serum samples demonstrating gross lipemia, gross hemolysis, or turbidity should not be used with this test.
4. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

**QUALITY CONTROL**

Good laboratory practice requires that quality control specimens be run with each calibration curve to verify assay performance. To assure proper performance, a statistically significant number of controls should be assayed to establish mean values and acceptable ranges. Controls containing sodium azide should not be used.

**SYMBOLS**

|   |                                    |
|---|------------------------------------|
|  | BATCH CODE                         |
|  | USE BY                             |
|  | MANUFACTURER                       |
|  | CONTAINS SUFFICIENT FOR <n> TESTS  |
|  | IN VITRO DIAGNOSTIC MEDICAL DEVICE |
|  | TEMPERATURE LIMITATION             |
|  | CATALOGUE NUMBER                   |
|  | CONSULT INSTRUCTIONS FOR USE       |

**REFERENCES**

1. Gerich JE: Hormonal control of homeostasis. IN Galloway JA, Potvin JH, Shuman CR: Diabetes Mellitus, ninth edition. Eli Lilly Co, Indianapolis, 1988 pp 46-63
2. Rasmussen H, Zawalick KC. Ganesan S, Calle R, Zawalich WS: Physiology and pathophysiology of insulin secretion. Diab Care 13:655-666, 1990
3. Gammeltoft S: Insulin receptors: binding kinetics and structure-function relationship of insulin. Physiol Rev 64:1321-1378, 1984
4. Rosen OM: After insulin binds, Science 237:1452-1458, 1987
5. Schwartz MW, Figlewicz D, Baskin DG, Woods SC, Porte D Jr. Insulin in the brain: A hormonal regulator of energy balance. Endocrin Rev 13:387-414, 1992

**For order and inquires, please contact**

AUTOBIO DIAGNOSTICS CO., LTD.  
ADD: No.87 Jingbei Yi Road, National Eco & Tech Development Area,  
Zhengzhou , China 450016  
Tel: +86-371-67985313 Fax: +86-371-67985804  
Web: [www.autobio.com.cn](http://www.autobio.com.cn)